

leCampion College

5th Form Laboratory

Skill: A/I /10

Name: Jada Dixon

Date: October 17, 2020

Topic: Volumetric Analysis

Aim: To determine the amount of sodium carbonate needed to standardize hydrochloric acid.

Title: Standardizing HCl using 0.1M Na₂CO₃ with methyl orange indicator

Apparatus and materials: burette, pipette, methyl orange, hydrochloric acid, sodium carbonate, conical flask,

Method:

1. The burette was conditioned and filled with HCl.
2. The pipette was conditioned with sodium carbonate
3. Using the pipette, 20 cm³ of 0.1M sodium carbonate was transferred to a conical flask.
4. Three drops of methyl orange were added to the conical flask.
5. Hydrochloric acid was added to the conical flask until one drop changed the colour from yellow to red.
6. The final volume of hydrochloric acid on the burette was recorded
7. The titration was repeated until the values were within 0.1cm³ of each other.

Results:

	Trial 1	Trial 2	Trial 3
Final Burette reading/cm ³	25.55	27.50	28.10
Initial Burette Reading/cm ³	3.25	6.55	7.25
Volume of HCl used /cm ³	22.30	20.95	20.85

Discussion

1. Write the chemical equation for the reaction [2]
 $2\text{HCl}(\text{aq}) + \text{Na}_2\text{CO}_3(\text{aq}) \rightarrow 2\text{NaCl}(\text{aq}) + \text{H}_2\text{O}(\text{l}) + \text{CO}_2(\text{g})$

2. Calculate the average volume of HCl

$$(20.95+20.85)/2$$

$$=20.90\text{cm}^3 \quad [1]$$

3. Calculate the moles of Na₂CO₃ in the pipetted volume
Concentration = moles/volume
Moles = concentration * volume
Moles = .1mol/dm³ * .02dm³
=.002 moles

[2]

4. Use the mole ratio to determine the moles of HCl in the average volume [1]

Mol ratio $\text{HCl}:\text{Na}_2\text{CO}_3 = 2:1$
X: .002

$X = 2 \times .002$

X=.004 moles HCl

5. Calculate the molar concentration of HCl
Concentration = #moles/volume
Concentration = $0.004\text{mol}/.0209\text{dm}^3$
Concentration = $0.191\text{mol}/\text{dm}^3$
[2]
6. Calculate the mass concentration of HCl
Mass concentration = molar concentration * molar mass
Mass concentration = $.191\text{mol}/\text{dm}^3 * 36.46\text{g}/\text{mol}$
Mass concentration = $6.96\text{g}/\text{dm}^3$

[2]

Conclusion: An average of 20.90cm^3 of HCl was needed to neutralize 20cm^3 of sodium carbonate of concentration 0.1M. There were .002 moles of sodium carbonate, and .004 moles of hydrochloric acid, which had a molar concentration of $.191\text{mol}/\text{dm}^3$ and a mass concentration of $6.96\text{g}/\text{dm}^3$.